

Agilent IDP Series Dry Scroll Pumps

Dry. High Performance.



High Performance. Compact. Dry.

High performance in a compact package – the Agilent IDP Series Dry Scroll Pumps provide affordable oil-free vacuum with easy system integration and they are suitable for a wide variety of applications.


The IDP Series employs an innovative hermetic design in which the motor and bearings are outside the vacuum space allowing full isolation of all pumped gases.

In addition, the IDP dry pump technology is environmentally friendly, as it eliminates the need for oil disposal or the handling of contaminated oil in your laboratory.

IEC-320 Power Connector and Switch




Hermetically Sealed, Fully Isolated Bearing



Scroll




Gas Ballast



Removable Polypropylene Exhaust Fitting



Understanding Scroll Technology



Gas enters scroll set



Gas is displaced and...



...compressed toward center hub



Gas exhausted at center hub

IDP Pumps create vacuum using a simple dual scroll mechanism in which one of the nested scrolls orbits about the other, creating moving zones of captured gas. Gas enters the scroll set at the perimeter and is displaced and compressed toward the center hub where it is exhausted.

The uncomplicated dual scroll design offers many benefits including lower noise and vibration levels, simple and infrequent maintenance requirements, plus the elimination of catastrophic failure modes.

High Pumping Speed. Low Base Pressure.

Agilent IDP Series Scroll Pumps, with high pumping speed and very low base pressure, provide all the advantages of scroll pump technology in a clean, compact, lightweight, and cost-effective package.



IDP-2 Dry Scroll Pump

- Pumping speed of 35 l/m (2.1 m³/hr)
- Achieves a base pressure of 750 milliTorr (1.0 mbar)
- Economical solution for the wide range of less demanding applications



IDP-3 Dry Scroll Pump

- Pumping speed of 60 l/m 7(3.6 m³/hr)
- Very low base pressure of less than 250 mTorr (0.3 mbar)
- Truly high-performance solution for demanding applications

Why Choose IDP Series?

IDP Series vs. Rotary Vane Pumps

- ✓ Oil-free technology of IDP Pumps eliminate the possibility of oil contamination in the vacuum system or of oil spills or leaks into the work environment.
- ✓ Maintenance of the IDP Pumps require only a simple, infrequent tip seal change as compared to oil checks, changes, and disposal.
- ✓ The IDP Pumps do not depend on the presence of sufficient oil to prevent seizing.

IDP Series vs. Membrane Diaphragm Pumps

- ✓ IDP Pumps are compact at 358 x 181 x 140 mm (14 x 7 x 6 in.), yet provide base pressure of less than 250 mTorr, almost 4 times lower than equivalently sized membrane/diaphragm pump.
- ✓ In turbo pump applications, the pumps' lower base pressure reduce power consumption and bearing temperature and increase the reliability of the system.
- ✓ IDP Pumps avoid catastrophic failure mechanisms. A diaphragm-pumped system may suffer sudden, rapid loss of pressure when a membrane ruptures.
- ✓ IDP Pumps produce lower noise and vibration levels than diaphragm pumps, creating a quiet, pleasant work environment, and much lower contribution to system noise and vibration.

Features and Benefits



Powerful – High Performance

IDP Dry Scroll Pumps deliver better vacuum performance than other pumps of similar size.

- High pumping speed provides rapid pumpdown and high gas throughput
- Very low base pressure ensures optimal turbo pump performance with increased system reliability



Oil-free – Affordable Dry Vacuum

Convenient, low maintenance pumps provide affordable dry vacuum for research and industrial applications alike.

- IDP Pumps eliminate hydrocarbon contamination in the vacuum system
- No oil leaks into the work environment



Innovative – Hermetic Design

IDP Pumps fully isolate the bearings and motor from the vacuum space and contain all pumped gases.

- Allows recovery of precious process gases
- Prevents leakage of toxic gases



Economical – Low Cost of Ownership

Robust, efficient scroll pump technology has a long service interval resulting in very low cost of ownership.

- Maintenance time and costs of oil topping, replacement and disposal are eliminated
- Simple tip seal or module exchange requires less than 30 minutes



System Friendly – Easy to Integrate

With a small footprint, low weight, and low power requirements, IDP Pumps are easy to accommodate in system design.

- IDP Pumps place little additional burden on system utilities and are suitable for use inside cabinet enclosures
- At only 21 pounds (9.5 kg), IDP Pumps are lighter than mechanical rotary vane and membrane pumps of similar pumping speed
- Low noise and vibration provide a quiet, pleasant work environment

Applications



Analytical Instrumentation

Monitoring, maintenance and environmental challenges of oil disposal make oil-sealed rotary vane pumps, at best, a nuisance in the laboratory environment. Without oil, high capacity IDP Pumps provide your instrument the required pumping speed and base pressure not available with other dry pump technologies.

- GC/MS
- Residual gas analyzers
- Load locks for mass spectrometers
- Electron microscopes
- Surface area analyzers
- EM Sample prep coaters



Semiconductor Fabrication

IDP Series hermetic design ensures against gas leakage into the environment, essential for evacuation and cycle purging of gas transfer lines and gas panels in wafer fabs. The pumps maintain ultra-pure gases by reducing the number of purge/vent cycles, so that on a single pumpdown the pumps low base pressure eliminates 99.9% of the volume gas.

- Gas panels
- Gas delivery systems
- Double wall containment vessels



High Energy Physics

Chamber evacuation and backing of turbo pumps on beam lines are just a couple of the applications that benefit from the clean, dry pump technology of the IDP Series. A compact, rugged, high performance design makes it a superior choice over membrane pumps, and the 55 dBA noise spec makes it virtually unnoticeable in the work environment.

- Accelerators
- Beam lines
- Synchrotrons



Industrial processing

IDP pumps deliver clean vacuum critical for many precision components – optics, lasers, electron devices – with a hermetically sealed design with all bearings outside the vacuum space. With low base pressure and high pumping speed, the IDP Series is the superior choice compared to membrane pumps.

- Laser evacuation
- Cryogenic dewars
- Glove boxes
- Gas recovery and recirculation
- Cryopump roughing
- Tube processing
- Vacuum chucks




General Research & Development

A cost effective solution for clean vacuum when fast pumpdown times are desired, the IDP Series outperforms other dry technologies with higher retained pumping speed both at high pressure and at base vacuum.

- Experimental chambers
- Sample prep systems
- Surface analysis systems
- Backing turbo pump on beam lines

Technical Specifications

		IDP-2	IDP-3
Peak pumping speed	60 Hz	35 l/m, 2.1 m ³ /hr, 1.2 cfm	60 l/m, 3.6 m ³ /hr, 2.1 cfm
	50 Hz	30 l/m, 1.8 m ³ /hr, 1.1 cfm	50 l/m, 3.0 m ³ /hr, 1.8 cfm
Ultimate pressure		7.5 x 10 ⁻¹ torr (1.0 mbar, 100 Pa)	2.5 x 10 ⁻¹ torr (3.3 x 10 ⁻¹ mbar, 33 Pa)
Maximum inlet pressure		1 atmosphere (1.0 bar, 101 kPa)	
Maximum outlet pressure		1.4 atmospheres (1.4 bar, 142 kPa)	
Inlet connection		NW16 KF flange	
Exhaust connection		Female 1/4 in. NPT (10 mm hose barb provided)	
Gas ballast connection		Female 1/8 in. NPT	
Ambient operating temperature		5 to 40°C (41 to 108°F)	
Storage temperature		-20 to 60 °C (-4 to 140°F)	
Motor rating		0.16 horsepower (0.12 kW)	
Supply power		1Ø – 100 VAC, 50-60 Hz : 115 VAC, 60 Hz : 220-240 VAC, 50-60 Hz	
Motor thermal protection		Automatic	
Rotation speed	60 Hz	3200 RPM	
	50 Hz	2600 RPM	
Cooling		Air-cooled	
Weight		9.5 kg (21 lbs)	Shipping-10.5 kg (23 lbs)
Restrictions		No corrosive, explosive, or particulate-forming gases	
Leak rate		<1 x 10 ⁻⁶ std-cc/sec helium	
Noise level (per ISO 11201)		55 dB(A)	
Vibration level at inlet (per ISO 10816-1)		1.5 mm/second	
Compliance		Conforms with CE, CSA, CSA/CUS, Semi S2-703, and RoHS	





Ordering Information

IDP Series Dry Vacuum Pumps		with Hour Meter	with Isolation Valve	with Hour Meter and Isolation Valve	Tip Seal Kit	Replacement Module
IDP-2, 1Ø, 220V, 50/60Hz	IDP2A01	IDP2A11	IDP2A21	IDP2A31	IDP2TS	IDP2
IDP-2, 1Ø, 115V, 50/60Hz	IDP2B01	IDP2B11	IDP2B21	IDP2B31		
IDP-2, 1Ø, 100V, 50/60Hz	IDP2C01	IDP2C11	IDP2C21	IDP2C31		
IDP-2, 24 VDC	IDP2D01	IDP2D11	IDP2D21	IDP2D31		
IDP-3, 1Ø, 220V, 50/60Hz	IDP3A01	IDP3A11	IDP3A21	IDP3A31	IDP3TS	IDP3
IDP-3, 1Ø, 115V, 50/60Hz	IDP3B01	IDP3B11	IDP3B21	IDP3B31		
IDP-3, 1Ø, 100V, 50/60Hz	IDP3C01	IDP3C11	IDP3C21	IDP3C31		
IDP-3, 24 VDC	IDP3D01	IDP3D11	IDP3D21	IDP3D31		

Power Cord Selection

Europe, 10A/220-230V, 2.5 Meter	656494220	Israel, 10A/230V, 2.5 Meter	656494230
Denmark, 10A/220-230V, 2.5 Meter	656494225	Japan, 12A/100V, 2.3 Meter	656494240
Switzerland, 10A/230V, 2.5 Meter	656494235	North America, 15A/125V, 2.0 Meter	656458203
UK/Ireland, 13A/230V, 2.5 Meter	656494250	North America, 10A/230V, 2.5 Meter	656494255
India, 10A/220-250V, 2.5 Meter	656494245		



Service & Support

United States and Canada

Agilent Technologies
121 Hartwell Avenue
Lexington, MA 02421 USA
Tel.: +1 781 861 7200
Toll-Free: +1 800 882 7426
Fax: +1 781 860 5437
vpl-customerservice@agilent.com

Benelux

Agilent Technologies Netherlands B.V.
Herculesweg 8
4338 PL Middelburg The Netherlands
Tel: +31 118 671570
Fax: +31 118 671569
Toll free: 00 800 234 234 00

China

Agilent Technologies (China) Co. Ltd
No.3, Wang Jing Bei Lu,
Chao Yang District,
Beijing, 100102 China
Tel.: +86 (10) 6439 7888
Fax: +86 (10) 6439 1318
Toll-Free: 800 820 8266
vpc-customerservice@agilent.com

France

Agilent Technologies France
7 avenue des Tropiques
Z.A. de Courtaboeuf – B.P. 12
91941 Les Ulis cedex France
Tel.: +33 (0) 1 69 86 38 84
Fax: +33 (0) 1 69 86 29 88
Toll free: 00 800 234 234 00
vpf.sales@agilent.com

Germany & Austria

Agilent Technologies
Alsfelder Strasse 6
Postfach 11 14 35
64289 Darmstadt Germany
Tel.: +49 (0) 6151 703 353
Fax: +49 (0) 6151 703 302
Toll free: 00 800 234 234 00

India

Agilent Technologies India Pvt. Ltd.
G01, Prime corporate Park,
230/231, Sahar Road,
Opp. Blue Dart Centre,
Andheri (East), Mumbai – 400 099 India
Tel: +91 22 30648287/8200
Fax: +91 22 30648250
Toll Free: 1800 113037
cag_india@agilent.com

Italy

Agilent Technologies Italia S.p.A.
via F.lli Varian 54
10040 Leini, (Torino) Italy
Tel.: +39 011 997 9111
Fax: +39 011 997 9350
Toll-Free: 00 800 234 234 00
vpt.sales@agilent.com
vpt-customerservice@agilent.com

Japan

Agilent Technologies Japan, Ltd.
8th Floor
Sumitomo Shibaura Building
4-16-36 Shibaura Minato-ku
Tokyo 108-0023 Japan
Tel.: +81 3 5232 1253
Toll-Free: 0120 655 040
Fax: +81 3 5232 1710
vpj-customerservice@agilent.com

Korea

Agilent Technologies Korea Ltd.
Shinsa 2nd Bldg. 2F
966-5 Daechi-dong
Kangnam-gu, Seoul Korea 135-280
Tel.: +82 2 3452 2455
Toll-Free: 080 222 2452
Fax: +82 2 3452 2451
vpk-customerservice@agilent.com

Mexico

Agilent Technologies
Concepcion Beistegui No 109
Col Del Valle
C.P. 03100 Mexico, D.F.
Tel.: +52 5 523 9465
Fax: +52 5 523 9472

Singapore

Agilent Technologies Singapore Pte Ltd
No.1 Yishun Avenue 7 Singapore 768923
Tel: +65 6215 8045
Fax: +65 6754 0574
Toll-Free: 1 800 2762622
vps-customerservice@agilent.com

South East Asia

Agilent Technologies Sales Sdn Bhd
Unit 201, Level 2 uptown 2,
2 Jalan SS21/37, Damansara Uptown
47400 Petaling Jaya, Selangor, Malaysia
Tel: +603 7712 6106
Fax: +603 6733 8121
Toll-Free: 1 800 880 805
vps-customerservice@agilent.com

Taiwan

Agilent Technologies Taiwan Limited
20 Kao-Shuang Rd.,
Pin-Chen City, 324
Taoyuan Hsien, Taiwan, R.O.C.
Tel. +886 34959281
Toll Free: 0800 051 342
vpw-customerservice@agilent.com

UK & Ireland

Agilent Technologies UK Ltd
6 Mead Road
Oxford Industrial Park
Yarnton, Oxford OX5 1QU UK
Tel.: +44 (0) 1865 291570
Fax: +44 (0) 1865 291571
Toll free: 00 800 234 234 00
vpt-customerservice@agilent.com

This information is subject to change without notice.

© Agilent Technologies, Inc., 2011
Published in USA, May, 2011



Agilent Technologies

Learn more:
www.agilent.com/chem/vacuum