




**Helium-calibrated test
leak with ISOlution
valve, leak rate between
 1 and $3 \cdot 10^{-7}$ mbar l/s,
aluminum reservoir,
DN 25 ISO-KF**

Helium Calibrated Test Leak with Isolation Valve, Leak Rate between 1 and 3 · 10⁻⁷ mbar l/s



Helium-calibrated test leak with ISolation valve, leak rate between 1 and 3 · 10⁻⁷ mbar l/s, aluminum reservoir, DN 25 ISO-KF

Dimensions



Order number	Helium-calibrated test leak with ISolation valve, leak rate between 1 and 3 · 10 ⁻⁷ mbar l/s, aluminum reservoir, DN 25 ISO-KF
Helium-calibrated test leak with ISolation valve, leak rate between 1 and 3 · 10 ⁻⁷ mbar l/s, aluminum reservoir, DN 25 ISO-KF	FV4710



Are you looking for an optimum vacuum solution?

Talk to us:
Pfeiffer Vacuum GmbH
Germany
T +49 6441 802-0

Or scan the barcode, to visit our web page:



<https://webportal.pfeiffer-vacuum.com/global/en/contact>



Errors and/or changes excepted. - 11/1/2024

Follow Us On Social Media
#pfeiffervacuum



www.pfeiffer-vacuum.com

