











Mounting kit for HiPace® 60/80

PFEIFFER



Illustration only

Mounting kit for HiPace 60/80, DN 63 ISO-K, including coated centering ring and claws $\,$

- For DN 63 ISO-K
- Suitable for HiPace 60/80 and SplitFlow 50
- Centering ring with multi-functional coating for anti-twist protection
- 4 x claws M8
- Note: Required for secure mounting of the turbopump

Technical Data	Mounting kit for HiPace 60/80, DN 63 ISO-K, including coated centering ring and claws
Number of lugs	(M8) 4 pcs.
Screw	M8
Connection flange	DN 63 ISO-K

	Mounting kit for HiPace 60/80, DN 63 ISO-K, including coated centering ring and claws
Mounting kit for HiPace 60/80, DN 63 ISO-K, including coated centering ring and claws	PM 016 510 -T







Are you looking for an optimum vacuum solution?

Talk to us: Pfeiffer Vacuum GmbH Germany T +49 6441 802-0

Or scan the barcode, to visit our web page:



https://webportal.pfeiffer-vacuum.com/global/en/contact

Follow Us On Social Media #pfeiffervacuum











www.pfeiffer-vacuum.com

