




Mounting kit for HiPace 60/80, DN 63 ISO-K, including coated centering ring, splinter and claw grips



PFEIFFER 

Mounting kit for HiPace® 60/80



Illustration only

Mounting kit for HiPace 60/80, DN 63 ISO-K, including coated centering ring, splinter and claw grips

- For DN 63 ISO-K
- Suitable for HiPace 60/80 and SplitFlow 50
- Centering ring with multi-functional coating for anti-twist protection and integrated splinter shield
- Splinter shield: to safeguard against foreign objects, mesh width 0.8 mm, molecular conductance for N₂ 194 l/s
- 4 x claws M8
- Note: Required for secure mounting of the turbopump

Technical Data	Mounting kit for HiPace 60/80, DN 63 ISO-K, including coated centering ring, splinter and claw grips
Number of lugs	4 pcs.
Molecular conductance for N ₂	194 l/s
Connection flange	DN 63 ISO-K

Order number	Mounting kit for HiPace 60/80, DN 63 ISO-K, including coated centering ring, splinter and claw grips
Mounting kit for HiPace 60/80, DN 63 ISO-K, including coated centering ring, splinter and claw grips	PM 016 511 -T



Errors and/or changes excepted. - 2/21/2025

Are you looking for an optimum vacuum solution?

Talk to us:

Pfeiffer Vacuum GmbH

Germany

T +49 6441 802-0

Or scan the barcode, to visit our web page:



<https://webportal.pfeiffer-vacuum.com/global/en/contact>

Follow Us On Social Media

#pfeiffervacuum



www.pfeiffer-vacuum.com

PFEIFFER 